

# RecombiLISA

## HBsAb ELISA

IVD REF E0711

- 96-well ELISA kit for the qualitative detection of hepatitis B surface antibody (HBsAb) in human serum or plasma
- For export only, not for re-sale in the USA
- Store at 2-8°C upon receipt

### INTENDED USE

The *RecombiLISA* HBsAb ELISA is a solid-phase enzyme-linked immunosorbent assay for the qualitative detection of hepatitis B surface antibody (HBsAb) in human serum or plasma. It is intended for professional use only as an aid in the diagnosis of infection with hepatitis B virus (HBV).

### SUMMARY AND EXPLANATION OF THE TEST

Hepatitis virus B is the most common cause of persistent viremia and the most important cause of chronic liver disease and hepatocellular carcinoma. Clinically apparent HBV infections may have been extant for several millennia. It is estimated that there are 300 million chronic carriers of HBV in the world. The carrier rates vary from as little as 0.3% (Western countries) to 20% (Asia, Africa)<sup>1</sup>.

HBV is a hepatotropic DNA virus. The core of the virus contains a DNA polymerase<sup>2</sup>, the core antigen (HBcAg)<sup>3</sup> and the e antigen (HBeAg)<sup>4</sup>. The core of HBV is enclosed in a coat that contains lipid, protein and carbohydrate and expresses an antigen termed hepatitis B surface antigen (HBsAg)<sup>3</sup>.

HBsAg is the first marker to appear in the blood in acute hepatitis B, being detected 1 week to 2 months after exposure and 2 weeks to 2 months before the onset of symptoms. Simultaneous with, or shortly after the disappearance of HBsAg, HBsAb is found in the blood. Its appearance heralds complete recovery, and its presence provides lifelong immunity<sup>1-7</sup>.

### TEST PRINCIPLE

The *RecombiLISA* HBsAb ELISA is a solid-phase enzyme-linked immunosorbent assay based on the principle of the double antigen sandwich technique for the detection of HBsAb in human serum or plasma.

The *RecombiLISA* HBsAb ELISA is composed of two key components:

- Solid microwells pre-coated with HBsAg
- Liquid conjugates composed of HBsAg conjugated with horseradish peroxidase (HRP-HBsAg Conjugate)

During the assay, the test specimen and HRP-HBsAg Conjugate are incubated simultaneously in the coated microwells. The HBsAb, if present in the specimen, binds to the HBsAg coated on the microwell surface as well as the antigen in the HRP-HBsAg Conjugate, forming a double antigen sandwich complex.

Unbound materials are then removed by washing. After addition of the TMB Substrate, the presence of the sandwich complex is shown by a blue color resulting from a reaction between the enzyme and substrate. The reaction is then quenched by addition of the Stop Solution and the absorbance value for each microwell is determined using a spectrophotometer at 450/620-690 nm.

### MATERIALS AND REAGENTS

#### Materials and reagents provided with the kit

Item	Description	Quantity	Catalog
1	HBsAg Coated Microwells	8 wells x 12 strips	E0711W
2	HBsAb Negative Control	1 mL	E0711N
3	HBsAb Positive Control	1 mL	E0711P
4	HRP-HBsAg Conjugate	6 mL	E0711H
5	Wash Buffer (30x Concentrate)	20 mL	WE3000
6	TMB Substrate A	6 mL	TME2000A
7	TMB Substrate B	6 mL	TME2000B
8	Stop Solution	6 mL	SE1000
9	ELISA Working Sheet	2	E0001ES
10	Product Insert	1	PI-E0711
Others	3 x Microplate Sealers and 1 x Resealable Bag		

#### Materials and reagents required but not provided in the kit

- Pipette capable of delivering 50 µL and 100 µL
- Microplate reader with a bandwidth of 10 nm or less and an optical density range of 0-3 OD or greater at 450 nm wavelength is acceptable
- Absorbent paper for blotting the microplate wells
- Distilled or de-ionized water
- Timer

### STORAGE AND STABILITY

All reagents except the concentrated wash buffer are ready to use as supplied. Store all components at 2-8°C. Do not freeze. Avoid strong light. Reseal the microwells after removing the desired number of wells. Place unused wells in the resealable bag provided and return to 2-8°C. Ensure that the reagents are brought to room temperature before opening. All the reagents are stable through the expiration date printed on the label if not opened.

### SPECIMEN COLLECTION AND PREPARATION

- Serum should be prepared from a whole blood specimen obtained by acceptable venipuncture technique.
- This kit is designed for use with serum specimen without additives only.
- If a specimen is not tested immediately, refrigerate at 2-8°C. If storage period greater than three days is anticipated, the specimen should be frozen (-20°C). Avoid repeated freezing-thawing of specimens. If a specimen is to be shipped, pack in compliance with federal regulation covering the transportation of etiologic agents.
- Specimens containing precipitants may give inconsistent test results. Clarify such specimens by centrifugation prior to assaying.
- Do not use serum specimens demonstrating gross lipemia, gross hemolysis or turbidity. Do not use specimens containing sodium azide.

### PREPARATION OF THE REAGENTS

- Bring all reagents, controls to room temperature (18-28°C).
- Preparation of working Wash Buffer.** Warm up the concentrated Wash Buffer to 37°C to dissolve the precipitant if it appears. Dilute concentrated Wash Buffer 30-fold with water as follows:

Plate	DI water	Wash buffer (30X)	Final volume
1 strip	58 mL	2.0 mL	60 mL
2 strips	116 mL	4.0 mL	120 mL
3 strips	174 mL	6.0 mL	180 mL
4 strips	232 mL	8.0 mL	240 mL

The diluted wash buffer can be stored at 2-8°C for up to 3 days.

- Mix each reagent before adding to the test wells.
- Determine the number of microwells needed and mark on the ELISA Working Sheet with the appropriate information. Positive and Negative Controls should be run in duplicate to ensure accuracy.











### ASSAY PROCEDURE

- Calculate the desired number of microwells. Remove the remaining microwells and place them with desiccant into the resealable plastic bag, seal and store at 2-8°C for later use.
- Add specimens according to the designation on the ELISA Working Sheet:
  - Blank Well:** Do not add any reagents.
  - Control Wells:** Add 50 µL of HBsAb Positive, Negative Control into the designated control wells, respectively.
  - Test Wells:** Add 50 µL of test specimen into each test well, respectively.
- Add 50 µL of HRP-HBsAg Conjugate into each well, but not blank well.
- Gently shake the wells for 20 seconds, and then cover the wells.
- Incubate the wells at 37°C for 30 minutes.
- Wash Step (Can be performed manually or with automated washing):  
**Manual washing:** Carefully remove the incubation mixture by disposing the solution into a waste container. Fill each well with 350 µL of diluted Wash Buffer and shake gently for 20-30 seconds. Discard the wash solution completely.

Repeat 4 more times. After completing the last wash step, tap the plate on absorbent paper to remove residual liquid.

**Automatic washing:** Automatic plate washer must be calibrated to ensure efficient washing. Aspirate incubation mixture from all wells completely. Fill each well with 350 µL of diluted wash buffer and soak for 20-30 seconds. Aspirate all wells completely. Repeat 4 more times.

- Add 50 µL of TMB Substrate A and 50 µL of TMB Substrate B into each well, including the Blank Well. Cover the plate with a sealer.
- Incubate at 37°C in dark for 10 minutes.
- Stop the reaction by adding 50 µL of Stop Solution to each well. Gently mix for 30 seconds. **It is important to make sure that all the blue color completely changes to a yellow color.**
- Set the microplate reader wavelength at 450 nm. Measure the absorbance (OD) of each well against the Blank Well within 15 minutes after adding Stop Solution. A filter of 620-690 nm can be used as a reference wavelength to optimize the assay result.

Flow chart of assay procedure			
1.	Secure strips in microwell frame		Number of strips
2.	Add controls or specimens		50 µL
3.	Add HRP-HBsAg Conjugate		50 µL
4.	Gently shake		20 seconds
5.	Incubate		37°C, 30 minutes
6.	Wash: manual or automatic		5 times
7.	Add TMB Substrate A and B, respectively		50 µL + 50 µL
8.	Incubate in dark		37°C, 10 minutes
9.	Add Stop Solution. Gently mix		50 µL 30 seconds
10.	Read result		450/620-690 nm within 15 minutes

### INTERPRETATION OF RESULTS

- Set up the cut-off value**  
The cut-off value =  $N \times 2.1$   
N: Mean OD of the negative control. Use  $N=0.05$  for calculation of the cut-off value if the mean OD is less than 0.05.
- Calculation of specimen OD ratio**  
Calculate an OD ratio for each specimen by dividing its OD value by the Cut-off Value as follows:  
$$\text{Specimen OD ratio} = \frac{\text{Specimen OD}}{\text{Cut-off Value}}$$
- Assay validation**  
The mean OD value of the HBsAb positive controls should be  $\geq 0.80$ .  
The mean OD value of the HBsAb negative controls should be  $\leq 0.10$ .  
  
Check the assay procedure including incubation time and temperature and repeat assay if above criteria is not met.
- Interpretation of the results**

	<u>Specimen OD ratio</u>
Negative	< 1.00
Positive	$\geq 1.00$

- The negative result indicates that there is no detectable HBsAb in the specimen.

- Specimens with OD ratio  $\geq 1.00$  are initially considered to be positive by the *RecombiLISA* HBsAb ELISA. They should be retested in duplicate before a final interpretation is made.
- Results within 10% of the cut-off value should be interpreted with caution (it is advisable to retest in duplicate the corresponding specimens when it is applicable).  
If after retesting the absorbance of one of the duplicates is equal or greater than the cut-off value, the initial result is repeatable and the specimen is considered to be positive with the *RecombiLISA* HBsAb ELISA, subject to the limitation of the procedure, described below.  
If after retesting of a specimen, the absorbance value of the 2 duplicates are less than the cut-off value, the initial result is non-repeatable and the specimen is considered to be negative with the *RecombiLISA* HBsAb ELISA.  
Non-repeatable reactions are often caused by:
  - Inadequate microwell washing
  - Contamination of negative specimens by serum or plasma with a high antibody titer
  - Contamination of the TMB Substrate by oxidizing agents (bleach, metal ions, etc.)
  - Contamination of the Stop Solution

#### WARNING AND PRECAUTIONS

##### For in Vitro Diagnostic Use

- This package insert must be read completely before performing the test. Failure to follow the insert gives inaccurate test results.
- Do not use expired kits.
- Bring all reagents to room temperature (18-28°C) before use.
- Do not use the components in any other type of test kit as a substitute for the components in this kit.
- Do not use serum derived from hemolyzed blood specimen for testing.
- Do not ingest the reagents. Avoid contact with eyes, skin and mouth. Wear protective clothing and disposable gloves while handling the kit reagents and clinical specimens. Wash hands thoroughly after performing the test.
- Do not smoke, drink, or eat in areas where specimens or kit reagents are being handled.
- Users of this test should follow the US CDC Universal Precautions for prevention of transmission of HIV, HBV and other blood-borne pathogens.
- Dispose of all specimens and materials used to perform the test as bio-hazardous waste.
- At the beginning of each incubation, and after adding Stop Solution, gently rock the microwells to ensure thorough mixing. Avoid the formation of air bubbles which results in inaccurate absorbance values. Avoid splashing liquid while rocking or shaking the wells.
- Do not allow the microwells to dry between the end of the washing operation and the reagent distribution.
- The enzyme-substrate reaction is very sensitive to metal ions. Thus, do not allow any metal element to come into contact with the conjugate or TMB Substrate.
- The TMB Substrate must be colorless. The appearance of color indicates that the reagent cannot be used and must be replaced. The TMB Substrate B must be stored in the dark.
- Use a new distribution tip for each specimen. Never use the specimen container to distribute conjugate and TMB Substrate.
- The wash procedure is critical. Wells must be aspirated completely before adding the Wash Buffer or liquid reagents. Automatic washers must be validated with the test kit prior to use. Insufficient washing will result in poor precision and falsely elevated absorbance values.**
- The microplate reader must be calibrated as per the manufacturer's instructions to ensure accurate determination of absorbance. A non-calibrated reader may lead to invalid test results.
- Avoid exposure to strong light during color development.

#### LIMITATION OF THE TEST


















- The Assay Procedure and the Assay Result Interpretation must be followed closely when testing the presence of HBsAb in serum or plasma from individual subjects. Failure to follow the procedure may give inaccurate results.

- The *RecombiLISA* HBsAb ELISA is limited to the qualitative detection of HBsAb in human serum or plasma. The intensity of the color does not have linear correlation with the antibody titer in the specimen.
- A negative result for an individual subject indicates absence of detectable HBsAb. However, a negative test result does not preclude the possibility of exposure to or infection with HBV.
- A negative result can occur if the quantity of HBsAb present in the specimen is below the detection limits of the assay, or the antibody that are detected are not present during the stage of disease in which a specimen is collected.
- Some specimens containing unusually high titer of heterophile antibodies or rheumatoid factor may affect expected results.
- Any use or interpretation of this test's results must also rely on other clinical findings and the professional judgment of health care providers.

#### REFERENCES

- Rubin. E. Acute and chronic viral hepatitis. Federation Proceedings. 38 (13): 2665-2673 (1979).
- Aach R.D., Grisham J.W., Paker C.W. Detection of Australia antigen by radioimmunoassay. Proc. Natl. Acad. Sci. USA 68:1056 (1971).
- Kim. C.Y. and Tikes. J.G. Purification and biophysical characterization of hepatitis antigen. J. Clin Invest. 52: 1176-1186 (1973).
- Ling C. M. Radioimmunoassay for hepatitis B virus. Manual of clinical immunology. Rose N. and Friedman H. editors.
- Caldwell. C.W., Barrel.J.T. Enzyme immunoassay for hepatitis B and its comparison to other methods. Clin Chem Acta. 81:305 (1977).
- Wolters G. et al., Enzyme linked immunosorbent assay for hepatitis B surface antigen. J. Infect. Dis. 136:311 (1977)
- Magnius L.O., Lindhoilil, A Lunchn, P and Carson. S.A. New antigen-antibody system. Clinical significance in long-term, carriers of hepatitis B surface antigen. J. Am. Med. Asso. 231:356-359 (1975).

#### Index of Symbols

	See instructions for use		Use by
	For in vitro diagnostic use only		Conjugate
	Catalog #		Coated microwells
	Lot number		Positive control
	Tests per kit		Negative control
	Do not reuse		TMB substrate
	Manufacturer		Wash buffer
	Date of manufacture		Stop solution
	Store between 2-8°C		

  
**CTK Biotech, Inc.**  
 13855 Stowe Drive,  
 Poway, CA 92064, USA  
 Tel: 858-457-8698  
 Fax: 858-535-1739  
 E-mail: info@ctkbiotech.com

PI-E0711 Rev. CP  
 Date released: 2018-11-27  
 English version