

**Hunan Kangxin Biotechnology Co., Ltd.
NT-proBNP Test kit
(Microfluidic Fluorescence Immunoluminescence Assay)
Performance evaluation comparative study data**

Chengdu VACURE Biotechnology Co., Ltd.
June 2022

CHENGDU VACURE Biotech CO.,LTD

ADD: 4th Floor, Building 7, Sanyi Innovation Center (Phase III), No. 733, Section 2, Furong Avenue, Wenjiang District, Chengdu, SiChuan Province, China.

info@vacurebiotech.com

www.vacurebiotech.com

+86 28 85155537

1 Purpose & Overview

Purpose: To evaluate and analyze the performance of fluorescence immunoassay for the determination of Amino-Terminal pro-brain natriuretic peptide (NT-proBNP) in human serum. The intra-assay precision CV was calculated with reference to EP15-A2. Refer to the EP9-A2 file to calculate the linear regression equation and correlation coefficient between Hunan Kangxin Biotechnology Co., Ltd. (hereafter Kangxin) and Roche diagnostics GmbH (hereafter Roche). In order to analyze whether the performance of Kangxin is good compared with domestic similar products.

2 Reagents&Instruments

(1) Instruments: Electrochemiluminescence automatic immunoassay analyzer cobas e 411

Reagents: Brain Natriuretic Titanium Precursor Detection Kit

(Electrochemiluminescence)

Manufacturer: Roche Diagnostics GmbH

(2) Instruments: Fluorescence Immunoassay Analyzer LYOFIA-I

Reagents: N-terminal brain natriuretic peptide precursor (NT-proBNP) Test Kit (Microfluidic Fluorescence Immunoluminescence Assay)

Manufacturer: Hunan Kangxin Biotechnology Co., Ltd.

Lot: 0216401

3 Test content

3. 1 Precision assessment

Test samples at two concentration levels in 2 detection systems, repeat the test 10 times, calculate the average of 10 test results ($\bar{X}$) and standard deviation (S), get the coefficient of variation (CV), the results should not be less than the value declared by the manufacturer.

Manufacturer(Kangxin) declared value: Kangxin (CV) : $\leq 10\%$;

3. 2 Comparison of system results

Refer to the method in EP9-A2 Method Comparison and Bias Assessment Using Patient Samples "Method Comparison and Bias Assessment Using Patient Samples" to measure samples on two systems respectively, each sample was measured once, and statistical analysis was performed on the detection data.

Kangxin's linearity is: 30-35000pg/mL, the methodological comparison of the linearity of 30-35000pg/mL and 30-1000pg/mL of the Kangxin's linearity was carried out respectively.

Taking the detection result of the comparison system as the X axis and the detection result of the test system as the Y axis, make a regression curve to obtain the regression formula and the correlation coefficient r.

4 Test Results

4. 1 Precision assessment

Repeat times (pg/mL)	Roche		Kangxin	
	341	1003	341	1003
1	345	1013	332	1008
2	349	987	342	991
3	346	997	347	983
4	347	1012	335	1007
5	345	985	340	1022
6	345	1001	325	1027
7	340	1017	302	1043
8	331	1015	299	1065
9	333	1006	293	1078
10	327	1001	304	1087
Average Value	340.8	1003.4	321.9	1031.1
Standard Deviation	7.69993	11.31567	20.33306	36.13386
CV	2.3%	1.1%	6.3%	3.5%

4. 2 Methodological comparison


Sample/Unit (pg/mL)	Roche	Kangxin
1	26	< 30
2	29	43
3	38	46
4	40	67
5	41	< 30
6	56	65
7	59	160
8	74	69
9	80	106
10	85	83
11	85	57
12	89	67
13	129	117
14	134	251
15	137	61
16	138	105
17	144	138

18	147	151
19	158	175
20	166	199
21	175	142
22	175	132
23	175	158
24	197	152
25	204	208
26	227	137
27	227	203
28	269	298
29	309	346
30	313	359
31	434	353
32	577	623
33	606	595
34	612	539
35	662	606
36	798	681
37	863	968
38	896	893
39	901	985
40	954	809
41	1056	1302
42	1284	1412
43	1541	1480
44	1988	2273
45	2141	2130
46	2238	2008
47	2496	2231
48	2503	2364
49	3014	3079
50	3219	3342
51	3553	3669
52	4053	4075
53	4161	4739
54	4379	4254
55	4413	4299
56	4451	4720


57	4561	5400
58	4600	4475
59	4601	4680
60	6081	5217
61	6816	6483
62	7560	6073
63	7682	7016
64	10367	8128
65	16121	18549
66	17837	15728
67	18631	23706
68	22038	19066
69	24222	28051
70	34301	30143

4. 2. 1 Data Analysis

Test Range: 30-35000 pg/mL: Take Roche's detection value as X, and Kangxin's detection value as Y, to make a regression curve.



Test Range: 30-1000 pg/mL: Take the Roche detection value as X, and Kangxin's detection value as Y, to make a regression curve.



5 Conclusion

In this study, in the detection range of 30-35000 pg/mL and 30-1000 pg/mL, the correlations of Kangxin detection system are $R^2 = 0.9697$, $R^2 = 0.9547$ respectively, It has a good correlation with the Roche Cobas automatic electrochemiluminescence detection system. Precision meets manufacturer's declared standards, compared with the domestic similar amino-terminal brain natriuretic peptide precursor detection kit, it meets the clinical requirements.

Reporter :

Reviewer :

Date :