OXA-48 K-SeT

www.corisbio.com
IFU-58R1/FR/04

Fabricant:

Coris BioConcept
CREALYS Science Park
Rue Guillaume Fouquet, 11
5032 GEMBLoux
BELGIUM
Tel.: +32(0)81.719.917
Fax: +32(0)81.719.919
info@corisbio.com
Produit en BELGIQUE

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XII. PROBLEMES TECHNIQUES / RECLAMATIONS

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XIII. REFERENCES BIBLIOGRAPHIQUES

- DW Wareham and MH Abdul Momin.** Rapid Detection of Carbapenemases in Enterobacteriaceae: Evaluation of the RESIST-3 O.K.N (OXA-48, KPC, NDM) Multiplexed Lateral Flow Assay. *J Clin Microbiol.* 2017 Feb 1. pii: JCM.02471-16
- E Rubio, Y Zboromyrska, C Pitar, I Campo, I Alejo-Cancho, A Fasanella, A Vergara, F Marco, J Vila.** Evaluation of a rapid immunochromatographic test for the detection of OXA-48 carbapenemase. *Diagn Microbiol Infect Dis.* 2017 Mar;87(3):266-267
- F Koroska, S Göttig, M Kaase, J Steinmann, S Gatermann, J Sommer, T Wille, G Plum, A Hamprecht.** Comparison of phenotypic tests and an immunochromatographic assay and development of a new algorithm for OXA-48-like detection. *J Clin Microbiol.* 2016 Dec 28. Pii: JCM.01929-16
- F Pasteran, L Denorme, I Ote, S Gomez, D De Belder, Y Glupczynski, P Bogaerts, B Ghiglione, P Power, P Mertens, A Corso.** Rapid identification of OXA-48 and OXA-163 subfamily in carbapenem resistant gram-negative bacilli with a novel immunochromatographic lateral flow assay. *J Clin Microbiol.* 2016 Aug; 54(11):2832-2836
- D. Meunier, A. Vickers, R. Pike, R.L. Hill, N. Woodford and K.L. Hopkins.** Evaluation of the K-SeT R.E.S.I.S.T. immunochromatographic assay for the rapid detection of KPC and OXA-48-like carbapenemases. *J Antimicrob Chemother.* 2016 Aug; 71 (8):2357-9
- Wareham DW, Shah R, Betts JW, Phee LM, Abdul Momin MH.** Evaluation of an immunochromatographic Lateral Flow Assay (OXA-48 K-SeT) for the Rapid Detection of OXA-48-like Carbapenemases in Enterobacteriaceae. *J Clin Microbiol.* 2016 Feb;54 (2):471-3
- Fernández J, Fleites A, Rodcio MR, Vazquez F.** Evaluation of OXA-48 K-Se T: an immunochromatographic assay for rapid detection of OXA-48-producing Enterobacteriaceae. *Diagn Microbiol Infect Dis.* 2016 May;85 (1):12-5
- Dortet L, Jousset A, Sainte-Rose V, Cuzon G, Naas T.** Prospective evaluation of the OXA-48 K-SeT assay, an immunochromatographic test for the rapid detection of OXA-48-type carbapenemases. *J Antimicrob Chemother.* 2016 Jul;71 (7):1834-40
- Glupczynski Y, Evrard S, Ote I, Mertens P, Huang TD, Leclipteux T, Bogaerts P.** Evaluation of two new commercial immunochromatographic assays for the rapid detection of OXA-48 and KPC carbapenemases from cultured bacteria. *J Antimicrob Chemother.* 2016 May;71(5):1217-22
- B.M. Willey, X. Trimi, R. Ioboni, D.A. Boyd, G. Ricci, D.N. Grohn, D. Terenzi, A. Mazzulli, L. Mataseje, M. Mulvey, P. Lo, T. Mazzulli, S.M. Poutanen.** The Coris BioConcept OXA48 K-SeT Immuno-Chromatographic Assay Detects OXA48-type Carbapenemases with High Sensitivity and Specificity. 26th European Congress of Clinical Microbiology and Infectious Diseases, Amsterdam April 9-12, 2016
- A. Abulaila, F. Erdem, Z. Aktas and O. Oncul** Evaluation Comparison of a novel OXA-48 K-Set test and blue-carba test in detection of carbapenemase-producing enterobacteriaceae with using PCR as reference method. 26th European Congress of Clinical Microbiology and Infectious Diseases, Amsterdam April 09 – 12, 2016
- O. Karatuna, M. Kaya, I. Akyar** Evaluation of the performance of OXA-48 K-SeT immunochromatographic test for rapid identification of OXA-48 carbapenemase-producing Enterobacteriaceae. 26th European Congress of Clinical Microbiology and Infectious Diseases, Amsterdam April 09 – 12, 2016
- C.S. Nodari, A. Barth, C. Magagnin, A. Zavascki, A. Gales, C. Carvalhaes** OXA-370 is rapidly detected from different culture media using OXA-48 K-SeT® immunochromatography. 26th European Congress of Clinical Microbiology and Infectious Diseases, Amsterdam April 09 – 12, 2016
- P. Bogaerts, S. Evrard, G. Cuzon, TD. Huang, T. Naas and Y. Glupczynski** Specificity of the OXA-48 immunochromatographic K-SeT for the detection of OXA-48 like in *Shewanella* spp. 26th European Congress of Clinical Microbiology and Infectious Diseases, Infectious Amsterdam April 09 – 12, 2016
- A. Sarria, R. Gomez-Gil, G. Ruiz, MP. Romero, J. Garcia-Rodríguez** Preliminary study of the OXA-48 card letitest method for the direct detection of OXA-48 carbapenemase in blood and plates culture. 26th European Congress of Clinical Microbiology and Infectious Diseases, Amsterdam April 09 – 12, 2016
- L. Dortet, A. Jousset, V. Sainte-Rose, G. Cuzon, and T. Naas** Prospective evaluation of the OXA-48 K-SeT for the detection of OXA-48-type carbapenemase producing Enterobacteriaceae. 26th European Congress of Clinical Microbiology and Infectious Diseases, Amsterdam April 09 – 12, 2016
- P. Bogaerts, S. Evrard, M. Dozen, TD. Huang and Y. Glupczynski** Impact of the isolation medium for the detection of OXA-48 and KPC-producing Gram negative bacteria by immunochromatographic assays. 26th European Congress of Clinical Microbiology and Infectious Diseases, Amsterdam April 09 – 12, 2016
- I. De Toro Peinado, M^oC.M. Gradolph, R. S. Rodriguez, M. V. Troya, M^oP. B. Ruiz, B. Palop** A rapid test for detection of OXA-48 carbapenemase. 26th European Congress of Clinical Microbiology and Infectious Diseases, Amsterdam April 09 – 12, 2016
- Bogaerts, P** Validation of a new lateral flow assay for the detection of OXA-48-like-producing Enterobacteriaceae. 25th ECCMID, Copenhagen, April, 26th, 2015, Abstract # 0261.

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¹ Newcombe, Robert G. "Two-Sided Confidence Intervals for the Single Proportion: Comparison of Seven Methods," *Statistics in Medicine*, 17, 857-872 (1998).